

Food & beverage

VT547

VT547

Low Temperature Freezer -60°C

The VT low temperature freezers creates the possibility to maintain temperatures as low as -60°C. Supreme stability, reliability, user-friendliness and ease of cleaning make these freezers an ideal solution for laboratories and hospitals. Nature R refrigerant, cyclopentane insulation combined with recyclable materials make the freezers extremely eco-friendly.

🏆 **Energy Efficient**

🏆 **Plug & Play**

🏆 **Porthole**



VT547 | Gallery

VT547 |

Features	Value
No. of baskets	3 no.
Lock	✓

VT547 |

Dimensions	Value
Dimension	932x1664x665
Weight	92,4/108

VT547 |

Cooling	Value
Refrigerant / amount	
Energy	v. 20 °C, 5,775 kWh/24h
Temperature range	-40 to -60 °C

