

# Food & beverage


VT147

## VT147

### Low Temperature Freezer –60°C

The VT low temperature freezers creates the possibility to maintain temperatures as low as -60°C. Supreme stability, reliability, user-friendliness and ease of cleaning make these freezers an ideal solution for laboratories and hospitals. Nature R refrigerant, cyclopentane insulation combined with recyclable materials make the freezers extremely eco-friendly.

 **Energy Efficient**

 **Plug & Play**

 **Porthole**



## VT147 | Gallery

## VT147 |

Features	Value
No. of baskets	1 no.
Lock	✓

## VT147 |

Dimensions	Value
Dimension	885x725x605
Weight	59,8/55,8

## VT147 |

Cooling	Value
Refrigerant / amount	Nature R
Energy	v. 20 °C, 3,707 kWh/24h
Frequenze	Hz
Temperature range	-40 to -60 °C



