

# Food & beverage

VT206

## VT206

### Low Temperature Freezer –45°C

The VT low temperature freezers creates the possibility to maintain temperatures as low as -45°C. Supreme stability, reliability, user-friendliness and ease of cleaning make these freezers an ideal solution for laboratories and hospitals. R290 refrigerant, cyclopentane insulation combined with recyclable materials make the freezers extremely eco-friendly.

 **Energy Efficient**

 **Plug & Play**

 **Porthole**



## VT206 | Gallery

## VT206 |

Dimensions	Value
Dimension	885x923x605
Weight	63,4/58,4

Cooling	Value
Refrigerant / amount	R290
Energy	v. 20 °C, 3,14 kWh/24h
Frequenze	50 Hz
Temperature range	-25 to -45 °C



