

Food & beverage

VT146

VT146

Low Temperature Freezer –45°C

The VT low temperature freezers creates the possibility to maintain temperatures as low as -45°C. Supreme stability, reliability, user-friendliness and ease of cleaning make these freezers an ideal solution for laboratories and hospitals. R404a refrigerant, cyclopentane insulation combined with recyclable materials make the freezers extremely eco-friendly.

 **Energy Efficient**

 **Plug & Play**

 **Porthole**



VT146 | Gallery

VT146 |

Dimensions	Value
Dimension	885x725x605
Weight	55,4/51,4

Cooling	Value
Refrigerant / amount	R290
Energy	v. 20 °C, 1,865 kWh/24h
Frequenze	50 Hz
Temperature range	-25 to +45 °C

